

- d. This code **MAY** contain zero or one [0..1] @**displayName**="Expiration Date" (CONF:4437-3401).

The value of the Expiration Date as a time stamp.

3. **SHALL** contain exactly one [1..1] **value** with @**xsi:type**="TS" (CONF:4437-3395).

Note: Expiration Date as a time stamp

Expiration Date as a time stamp.

- a. This value **SHALL** contain exactly one [1..1] @**value** (CONF:4437-3402).

Figure 7: Expiration Date

```
<observation classCode="OBS" moodCode="EVN">
  <templateId root="2.16.840.1.113883.10.20.22.4.309" extension="2019-06-21"/>
  <code code="C101670" displayName="Expiration Date"
    codeSystem="2.16.840.1.113883.3.26.1.1" codeSystemName="NCI Thesaurus"/>
  <value xsi:type="TS" value="20221015"/>
</observation>
```

2.7 Implantable Device Status Observation

[observation: identifier urn:hl7ii:2.16.840.1.113883.10.20.22.4.305:2019-06-21 (open)]

Table 15: Implantable Device Status Observation Contexts

Contained By:	Contains:
UDI Organizer (optional)	

This template is intended to be used in addition to the **Product Instance** template urn:oid:2.16.840.1.113883.10.20.22.4.37 to augment the parsed data from the a Unique Device Identifier (UDI). This template is used to exchange the status of the patient's implantable medical device. This status is only relevant to medical devices implanted in the patient's body.

Table 16: Implantable Device Status Observation Constraints Overview

XPath	Card.	Verb	Data Type	CONF#	Value
observation (identifier: urn:hl7ii:2.16.840.1.113883.10.20.22.4.305:2019-06-21)					
templateId	1..1	SHALL		4437-3502	
@root	1..1	SHALL		4437-3505	2.16.840.1.113883.10.20.22.4.305
@extension	1..1	SHALL		4437-3506	2019-06-21
code	1..1	SHALL		4437-3503	
@code	1..1	SHALL		4437-3507	urn:oid:2.16.840.1.113883.3.26.1.1 (NCI Thesaurus (NCIt)) = C160939
@codeSystem	1..1	SHALL		4437-3508	2.16.840.1.113883.3.26.1.1
@codeSystemName	0..1	MAY		4437-3509	NCI Thesaurus
@displayName	0..1	MAY		4437-3510	Implantable Device Status
value	1..1	SHALL	CD	4437-3504	urn:oid:2.16.840.1.113762.1.4.1021.48 (Implantable Device Status)
@code	1..1	SHALL		4437-3511	
@codeSystem	1..1	SHALL		4437-3512	

1. **SHALL** contain exactly one [1..1] **templateId** (CONF:4437-3502).
 - a. This templateId **SHALL** contain exactly one [1..1] **@root**="2.16.840.1.113883.10.20.22.4.305" (CONF:4437-3505).
Note: template oid

Revision from "2018-05-01".

- b. This templateId **SHALL** contain exactly one [1..1] **@extension**="2019-06-21" (CONF:4437-3506).
2. **SHALL** contain exactly one [1..1] **code** (CONF:4437-3503).
Note: Code for "Implantable Device Status"

Code used to specify the **Implantable Device Status**.

- a. This code **SHALL** contain exactly one [1..1] **@code**="C160939" Implantable Device Status (CodeSystem: NCI Thesaurus (NCIt))
urn:oid:2.16.840.1.113883.3.26.1.1 **STATIC** (CONF:4437-3507).
- b. This code **SHALL** contain exactly one [1..1] **@codeSystem**="2.16.840.1.113883.3.26.1.1" (CONF:4437-3508).

- c. This code **MAY** contain zero or one [0..1] @codeSystemName="NCI Thesaurus" (CONF:4437-3509).

Implantable Device Status

- d. This code **MAY** contain zero or one [0..1] @displayName="Implantable Device Status" (CONF:4437-3510).
- 3. **SHALL** contain exactly one [1..1] value with @xsi:type="CD", where the code **SHALL** be selected from ValueSet Implantable Device Status
urn:oid:2.16.840.1.113762.1.4.1021.48 **STATIC** 2019-06-21 (CONF:4437-3504).
 - a. This value **SHALL** contain exactly one [1..1] @code (CONF:4437-3511).
 - b. This value **SHALL** contain exactly one [1..1] @codeSystem (CONF:4437-3512).

Figure 8: Implantable Device Status

```
<observation classCode="OBS" moodCode="EVN">
  <templateId root="2.16.840.1.113883.10.20.22.4.305" extension="2019-06-21"/>
  <code code="C160939" displayName="Implantable Device Status"
    codeSystem="2.16.840.1.113883.3.26.1.1" codeSystemName="NCI Thesaurus"/>
  <value xsi:type="CD" code="C160942" codeSystem="2.16.840.1.113883.3.26.1.1"
    displayName="Reduced Function" codeSystemName="NCI Thesaurus"
    sdct:valueSet="2.16.840.1.113762.1.4.1021.48">
    <originalText>
      <reference value="#ImplantableDeviceStatus_1"/>
    </originalText>
  </value>
</observation>
```

2.8 Latex Safety Observation

[observation: identifier urn:hl7ii:2.16.840.1.113883.10.20.22.4.314:2019-06-21 (open)]

Table 17: Latex Safety Observation Contexts

Contained By:	Contains:
UDI Organizer (optional)	

This template is intended to be used in addition to the **Product Instance** template urn:oid:2.16.840.1.113883.10.20.22.4.37 to exchange the **Latex Safety Status** of the patient's medical device. The UDI-DI of the medical device may be used to retrieve the Latex Safety status in accessGUDID, which should be considered the source of truth.