3. **SHALL** contain exactly one [1..1] **value** with @xsi:type="ED" (CONF:4437-3449). Note: Brand Name as a text value

Figure 4: Company Name

## 2.4 Device Identifier Observation

[observation: identifier urn:hl7ii:2.16.840.1.113883.10.20.22.4.304:2019-06-21 (open)]

Table 9: Device Identifier Observation Contexts

Contained By:	Contains:
<u>UDI Organizer</u> (required)	

This template is intended to be used in addition to the **Product Instance** template urn:oid:2.16.840.1.113883.10.20.22.4.37 to exchange the **Device Identifier** (also known as the "Primary DI Number") for a medical device marketed in the US. The device identifier is parsed from the UDI value.

The Device Identifier number can be used as a key to look-up device identification information in the publicly available version of the US FDA Global UDI Database (GUDID) - the AccessGUDID: https://accessgudid.nlm.nih.gov/resources/home

Table 10: Device Identifier Observation Constraints Overview

XPath	Card.	Verb	Data Type	CONF#	Value
observation (identifier: urn:hl7ii:	2.16.840.1.	113883.10.2	0.22.4.30	4:2019-06-2	1)
templateId	11	SHALL		4437- 3421	
@root	11	SHALL		4437- 3424	2.16.840.1.113883.10.20.22.4 .304
@extension	11	SHALL		4437- 3425	2019-06-21
code	11	SHALL		4437- 3422	
@code	11	SHALL		4437- 3426	urn:oid:2.16.840.1.113883.3.2 6.1.1 (NCI Thesaurus (NCIt)) = C101722
@codeSystem	11	SHALL		4437- 3427	2.16.840.1.113883.3.26.1.1
@codeSystemName	01	MAY		4437- 3428	NCI Thesaurus
@displayName	01	MAY		4437- 3429	Primary DI Number
value	11	SHALL	II	4437- 3423	
@root	11	SHALL		<u>4437-</u> <u>3430</u>	
@extension	11	SHALL		<u>4437-</u> <u>3431</u>	
@displayable	11	SHALL		<u>4437-</u> <u>3540</u>	true

- 1. **SHALL** contain exactly one [1..1] **templateId** (CONF:4437-3421).
  - a. This templateId **SHALL** contain exactly one [1..1] @root="2.16.840.1.113883.10.20.22.4.304" (CONF:4437-3424). Note: template oid

Revision from "2018-05-01"

- b. This templateId **SHALL** contain exactly one [1..1] **@extension=**"2019-06-21" (CONF:4437-3425).
- 2. **SHALL** contain exactly one [1..1] **code** (CONF:4437-3422). Note: Primary DI Number Code from NCIt

Code used to specify the device identifier, one of the parsed UDI components known as the Primary DI Number (UDI-DI).

a. This code **SHALL** contain exactly one [1..1] @code="C101722" Primary DI Number (CodeSystem: NCI Thesaurus (NCIt) urn:oid:2.16.840.1.113883.3.26.1.1 **STATIC**) (CONF:4437-3426).

- b. This code **SHALL** contain exactly one [1..1] @codeSystem="2.16.840.1.113883.3.26.1.1" (CONF:4437-3427).
- c. This code MAY contain zero or one [0..1] @codeSystemName="NCI Thesaurus" (CONF:4437-3428).
- d. This code **MAY** contain zero or one [0..1] **@displayName**="Primary DI Number" (CONF:4437-3429).

The value of the Device Identifier assigned using one of the FDA-accredited, UDI issuing agency: GS1, HIBCC, or ICBBA.

3. **SHALL** contain exactly one [1..1] **value** with @xsi:type="II" (CONF:4437-3423). Note: This value is assigned by the Manufacturer by using an FDA-accredited UDI-issuing

The "root" specifies the assigning authority, i.e, the UDI issuing agency. The allowable OIDs include:

- 1.3.160 (GS1)
- 2.16.840.1.113883.6.18 (ICCBBA ISBT-128)
- 2.16.840.1.113883.6.40 (HIBCC)
  - a. This value **SHALL** contain exactly one [1..1] **@root** (CONF:4437-3430).

The value of the Device Identifier.

b. This value **SHALL** contain exactly one [1..1] **@extension** (CONF:4437-3431).

The Device Identifier can be used to look up device information entered by manufacturers in the FDA GUDID

c. This value **SHALL** contain exactly one [1..1] **@displayable**="true" (CONF:4437-3540).

## Figure 5: Device Identifier

```
<observation classCode="OBS" moodCode="EVN">
   <templateId root="2.16.840.1.113883.10.20.22.4.304" extension="2019-06-21"/>
   <code code="C101722" displayName="Primary DI Number"</pre>
        \verb|codeSystem="2.16.840.1.113883.3.26.1.1"| | \verb|codeSystemName="NCI Thesaurus"| > \\
        <!-- GS1 Device Identifier -->
   <value xsi:type="II" root="1.3.160" extension="00999998989" displayable="true"</pre>
                  assigningAuthorityName="GS1"/>
</observation>
```