

- d. This code **MAY** contain zero or one [0..1] @**displayName**="Manufacturing Date" (CONF:4437-3467).

The value of the Manufacturing Date as a time stamp.

3. **SHALL** contain exactly one [1..1] **value** with @**xsi:type**="TS" (CONF:4437-3461).

Note: Manufacturing Date as a time stamp

Manufacturing Date as a time stamp.

- a. This value **SHALL** contain exactly one [1..1] @**value** (CONF:4437-3468).

Figure 11: Manufacturing Date

```
<observation classCode="OBS" moodCode="EVN">
  <templateId root="2.16.840.1.113883.10.20.22.4.316" extension="2019-06-21"/>
  <code code="C101669" displayName="Manufacturing Date"
    codeSystem="2.16.840.1.113883.3.26.1.1" codeSystemName="NCI Thesaurus"/>
  <value xsi:type="TS" value="20181015"/>
</observation>
```

2.11 Model Number Observation

[observation: identifier urn:hl7ii:2.16.840.1.113883.10.20.22.4.317:2019-06-21 (open)]

Table 22: Model Number Observation Contexts

Contained By:	Contains:
UDI Organizer (optional)	

This template is intended to be used in addition to the **Product Instance** template urn:oid:2.16.840.1.113883.10.20.22.4.37 to exchange the **Model Number** associated with the device. The UDI-DI of the medical device may be used to retrieve the **Model Number** in accessGUDID, which should be considered the source of truth.

Table 23: Model Number Observation Constraints Overview

XPath	Card.	Verb	Data Type	CONF#	Value
observation (identifier: urn:hl7ii:2.16.840.1.113883.10.20.22.4.317:2019-06-21)					
templateId	1..1	SHALL		4437-3412	
@root	1..1	SHALL		4437-3414	2.16.840.1.113883.10.20.22.4.317
@extension	1..1	SHALL		4437-3415	2019-06-21
code	1..1	SHALL		4437-3413	
@code	1..1	SHALL		4437-3416	urn:oid:2.16.840.1.113883.3.26.1.1 (NCI Thesaurus (NCIt)) = C99285
@codeSystem	1..1	SHALL		4437-3417	2.16.840.1.113883.3.26.1.1
@codeSystemName	0..1	MAY		4437-3418	NCI Thesaurus
@displayName	0..1	MAY		4437-3419	Model Number
value	1..1	SHALL	ED	4437-3420	

1. **SHALL** contain exactly one [1..1] **templateId** (CONF:4437-3412).

- a. This templateId **SHALL** contain exactly one [1..1]
@root="2.16.840.1.113883.10.20.22.4.317" (CONF:4437-3414).
 Note: template oid

Revision from "2018-05-01".

- b. This templateId **SHALL** contain exactly one [1..1] **@extension="2019-06-21"**
 (CONF:4437-3415).

2. **SHALL** contain exactly one [1..1] **code** (CONF:4437-3413).

Note: Code for "Device Model Number"

Code used to specify the **Model Number**.

- a. This code **SHALL** contain exactly one [1..1] **@code="C99285" Model Number**
 (CodeSystem: NCI Thesaurus (NCIt) urn:oid:2.16.840.1.113883.3.26.1.1
STATIC) (CONF:4437-3416).
- b. This code **SHALL** contain exactly one [1..1]
@codeSystem="2.16.840.1.113883.3.26.1.1" (CONF:4437-3417).
- c. This code **MAY** contain zero or one [0..1] **@codeSystemName="NCI Thesaurus"**
 (CONF:4437-3418).
- d. This code **MAY** contain zero or one [0..1] **@displayName="Model Number"**
 (CONF:4437-3419).

The Model Number can point to its corresponding narrative using the approach defined in CDA Release 2, section 4.3.5.1.

3. **SHALL** contain exactly one [1..1] **value** with @xsi:type="ED" (CONF:4437-3420).

Figure 12: Model Number

```
<observation classCode="OBS" moodCode="EVN">
  <templateId root="2.16.840.1.113883.10.20.22.4.317" extension="2019-06-21"/>
  <code code="C99285" displayName="Model Number"
    codeSystem="2.16.840.1.113883.3.26.1.1" codeSystemName="NCI Thesaurus"/>
  <value xsi:type="ED">
    <reference value="#ModelNumber_1"/>
  </value>
</observation>
```

2.12 MRI Safety Observation

[observation: identifier urn:hl7ii:2.16.840.1.113883.10.20.22.4.318:2019-06-21 (open)]

Table 24: MRI Safety Observation Contexts

Contained By:	Contains:
UDI Organizer (optional)	

This template is intended to be used in addition to the **Product Instance** template urn:oid:2.16.840.1.113883.10.20.22.4.37 to exchange the MRI Safety Status of the patient's medical device. The UDI-DI of the medical device may be used to retrieve the **MRI Safety Status** in accessGUDID, which should be considered the source of truth.